



STANDARD SPECIFICATION

Models	Chlorine and Sulfur Analyzer/Total Organic Halogen Analyzer Model TOX-300
Analysis methods	Oxidation decomposition / Coulometric titration
Oxidation decomposition	Quartz tube combustion method
Samples	Liquids, solids and gaseous (Use the gas injector Model GI-210.) samples
Sample Insertion	Automatic insertion by the sample boat (Use ABC unit.)
Pyrolytic furnace	Horizontal furnace: Up to 1100 degrees Celsius
Detection method	Oxidation-reduction potential (potentiometric detection by potential)
Detection electrode	Chlorine measurement: Silver electrode Sulfur measurement: Platinum electrode
Titration control	Automatic control of electrolytic current
Measurement range	Refer to Table Measurement Ranges".

Measurement Ranges

Absolute volume	Solid samples	Liquid samples	Gas samples
Chlorine: 0.05µg/ml to 50 µg	2µg/g (30mg)	0.5µg/ml (100µl)	5mg/m ³ (10ml)
Sulfur: 0.05µg/ml to 50 µg	2µg/g (30mg)	0.5µg/ml (100µl)	5mg/m ³ (10ml)

(Usual samples: Samples which are combustible under 1100 degrees Celsius and have no interference.)

Sample volume	Liquid sample: 100 µl or less, Solid sample: 30mg or less
Measurement Time	Within 10 minutes/ One measurement (At 2µg measurement)
Temperature and humidity	15 to 35 degrees Celsius, 80%RH or less (No condensation)
Dimensions	TOX-300 main unit: 550(W) × 360 (D) × 437 (H) mm ABC-210 unit: 440(W) × 250(D) × 180(H) mm
Weight	TOX-300 main unit: 36kg ABC-210 unit: 11kg
Computer	Windows® English version pre-install
Printer	Windows®-compatible printer
Gas	Oxygen (purity : 99.7 % or more) Argon (purity : 99.98 % or more)
Power	TOX-300 main unit: AC100/115/230/240V, 50/60Hz, 980VA ABC-210 unit: AC100/115/230/240V, 50/60Hz, 40VA



Note: Follow instructions in manuals to correctly install, connect and operate the instruments. Contents of catalogues are subject to change without prior notice when improvements are made in performance. The actual color of the goods may appear different from color printed. All screen images are simulated. *Company and product names contained herein are the trademarks or registered trademarks of the company concerned.

Your agent for Switzerland

Analytical instruments and services



Le Villaret
Rte. de la Robellaz 6
CH - 1432 Belmont s/Yverdon
Switzerland

Tel: +41 (0) 24 435 22 02
Fax: +41 (0) 24 435 22 03
Mobile: +41 (0) 79 435 96 26
Email: tecdev@swissonline.ch
Visit us: www.tecdevexpress.com

Nittoseiko Analytech Co.,Ltd.

7-10-1 Chuo-rinkan, Yamato, Kanagawa 242-0007, JAPAN
Tel: +81(0)46-278-0056
URL: <https://www.mccat.co.jp/global>
CAT No.05080520041E

Total Organic Halogen Analyzer Model TOX-300

*Cost effective solution for environmental analysis
Matrix independent combustion program without soot
Improve lower measuring limit*




Nittoseiko Analytech Co.,Ltd.

Advanced database software and the reputed coulometry methods allow measurement of chlorine and sulfur content in solid, liquid and gaseous samples down to the ppm level – to be fully utilized at various plants.
(Conforms to ASTM D3120,3246,3961,5808,JIS K2541)

Feature

1 Matrix Independent, unique single program for all samples.

Secure Combustion Program (SCP) enables ideal pyrolysis of substances in sample.



2 High-Speed Analysis Using the Boat Controller

Measurement is performed automatically in 6 to 12 minutes simply by placing a sample boat and starting measurement. Since the sample boat inlet box is provided with a cooling function (electronic cooling), the boat can be cooled in a short period of time for faster repetitive.




3 Application-Oriented System Configuration

The use of various options according to the particular sample make it possible to construct a system that precisely matches the application.

Basic Configuration

Electric furnace & Detector (horizontal type)
Control unit (PC and printer)

Options (Injector, Sampler, Other)




TOX-300 uses the same ABC and autosamplers as NSX-2100H. Customers can upgrade to high performance NSX-2100H by purchasing HF-210 and MCD-210.

Measuring Principle



Chlorine Analysis

Standard Set with Chlorine Titration Cell Unit
Samples are burned in an Argon/Oxygen atmosphere. The resulting chloride is lead into a titration cell where it is automatically titrated by silver ions generated coulometrically. The amount of chlorine is then calculated from the quantity of electricity required for the titration.



Sulfur Analysis

Standard Set with Sulfur Titration Cell Unit
Samples are burned in an Argon/Oxygen atmosphere. The resulting sulfur dioxide is lead into a titration cell where it is automatically titrated by triiodide ions generated coulometrically. The amount of sulfur is then calculated from the quantity of electricity required for the titration.



Application

Chlorine Analysis

Sample Name	Sample Size(mg)	No.of Measurements	Measured Value(ppm)	RSD(%)
Lubricating Oil	30	3	34	4
Kerosene	30µℓ	3	31	3
Crude Oil	30	3	4.0	3
Rubber	20	3	210	4
Carbon Fiber	20	3	193	3
Epoxy Resin	10	3	890	5
Polycarbonate Resin	30	3	7.9	3
Photoresist	30	3	7.2	2
Aluminum Foil	20	3	5.5	6
Waste Oil	10	3	390	3
Cement	10	3	460	3

Sulfur Analysis

Sample Name	Sample Size(mg)	No.of Measurements	Measured Value(ppm)	RSD(%)
Lubricating Oil	3	3	0.76%	3
Heavy Oil	10	3	1.2%	3
Engine Oil	5	3	0.41%	4
Polypropylene	30	3	2.0	5
Surface Activator	5	3	883	0.5
Rubber	15	3	0.55%	1
Manganese Oxide	30	3	37	7
Zinc Oxide	30	3	14	0.7
Lanolin	30	3	39	5
Fish Oil	20	3	180	4
Soil	5	3	369	6

Mode AOX

Total Organic Halogen Measurement (AOX)

Easy, rapid and accurate Total Organic Halogen measurement in environmental and industrial waste water by Mitsubishi adsorption and coulometry technique based on ISO9562, EPA9020 etc.

Measuring principle

1. Adsorption and washing

For adsorbing AOX, flow sample water, 10 – 200ml, through activated carbon with constant rate. After adsorption, wash the carbon by nitrate solution for eliminating inorganic halogens.

2. Combustion

The activated carbon on a sample boat is inserted into a heater furnace by Auto Boat Controller (ABC). AOX is converted to hydrogen halide by oxidative combustion.

3. Titration

Hydrogen halide is introduced into a titration cell for silver titration. Halogen content is calculated by electrical quantities based on Faraday's law.

Adsorption module TXA-04

Optional unit for automatic AOX adsorption and inorganic halogen washing.



Syringe	Disposable syringe 20ml
Flow method	Automatic syringe system
Maximum Flow Volume	
Syringe for flowing samples:	300ml (10 × n ml (n=1 to 30))
Syringe for washing:	50ml (5 × n ml (n=1 to 10))
Syringe quantity	5 pcs (for flowing samples or washing inorganic halide)
Column	3mm φ × 40mm, two-stage glass column
Power	AC 100V to 240V, 50/60Hz, 50VA
Dimensions / Weight	Approx. 480(W) × 270(D) × 530(H) mm / 13 kg

Option

	<table border="1"> <thead> <tr><th>MODEL</th><td>ASC-240S Solid Sample Changer</td></tr> <tr><th>Sample</th><td>Solid, Liquid (manual)</td></tr> <tr><th>Amount of sample</th><td>Solid 150mg Liquid 100µl</td></tr> <tr><th>Boat, number of sample</th><td>Ceramic, 40 pos.</td></tr> <tr><th>Sample injection</th><td>Auto boat control</td></tr> <tr><th>Boat cooling</th><td>Peltier</td></tr> <tr><th>Power</th><td>AC 100V to 240V, 50/60Hz, 80VA</td></tr> <tr><th>Dimensions / Weight</th><td>480(W) × 460(D) × 520(H) mm / 31 kg</td></tr> </thead> </table>	MODEL	ASC-240S Solid Sample Changer	Sample	Solid, Liquid (manual)	Amount of sample	Solid 150mg Liquid 100µl	Boat, number of sample	Ceramic, 40 pos.	Sample injection	Auto boat control	Boat cooling	Peltier	Power	AC 100V to 240V, 50/60Hz, 80VA	Dimensions / Weight	480(W) × 460(D) × 520(H) mm / 31 kg		<table border="1"> <thead> <tr><th>MODEL</th><td>ABC-210 Auto Boat Controller</td></tr> <tr><th>Sample</th><td>Solid, Liquid</td></tr> <tr><th>Amount of sample</th><td>Solid 150mg Liquid 100µl</td></tr> <tr><th>Boat</th><td>Quartz, disposable ceramic</td></tr> <tr><th>Boat cooling</th><td>Peltier</td></tr> <tr><th>Power</th><td>AC 100V to 240V, 50/60Hz, 40VA</td></tr> <tr><th>Dimensions / Weight</th><td>445(W) × 250(D) × 180(H) mm / 9 kg</td></tr> </thead> </table>	MODEL	ABC-210 Auto Boat Controller	Sample	Solid, Liquid	Amount of sample	Solid 150mg Liquid 100µl	Boat	Quartz, disposable ceramic	Boat cooling	Peltier	Power	AC 100V to 240V, 50/60Hz, 40VA	Dimensions / Weight	445(W) × 250(D) × 180(H) mm / 9 kg
MODEL	ASC-240S Solid Sample Changer																																
Sample	Solid, Liquid (manual)																																
Amount of sample	Solid 150mg Liquid 100µl																																
Boat, number of sample	Ceramic, 40 pos.																																
Sample injection	Auto boat control																																
Boat cooling	Peltier																																
Power	AC 100V to 240V, 50/60Hz, 80VA																																
Dimensions / Weight	480(W) × 460(D) × 520(H) mm / 31 kg																																
MODEL	ABC-210 Auto Boat Controller																																
Sample	Solid, Liquid																																
Amount of sample	Solid 150mg Liquid 100µl																																
Boat	Quartz, disposable ceramic																																
Boat cooling	Peltier																																
Power	AC 100V to 240V, 50/60Hz, 40VA																																
Dimensions / Weight	445(W) × 250(D) × 180(H) mm / 9 kg																																
	<table border="1"> <thead> <tr><th>MODEL</th><td>ASC-250L Liquid sample changer</td></tr> <tr><th>Sample</th><td>Liquid (non-aqueous, aqueous)</td></tr> <tr><th>Injection</th><td>Max 150µl (depend on sample)</td></tr> <tr><th>Injection speed</th><td>0.4 - 50µl/sec (depend on sample)</td></tr> <tr><th>Number</th><td>50pos in each 2, 4, 6ml vial tray</td></tr> <tr><th>Power</th><td>AC 100V to 240V, 50/60Hz, 180VA</td></tr> <tr><th>Dimensions / Weight</th><td>460(W) × 320(D) × 470(H) mm / 16 kg</td></tr> </thead> </table>	MODEL	ASC-250L Liquid sample changer	Sample	Liquid (non-aqueous, aqueous)	Injection	Max 150µl (depend on sample)	Injection speed	0.4 - 50µl/sec (depend on sample)	Number	50pos in each 2, 4, 6ml vial tray	Power	AC 100V to 240V, 50/60Hz, 180VA	Dimensions / Weight	460(W) × 320(D) × 470(H) mm / 16 kg		<table border="1"> <thead> <tr><th>MODEL</th><td>GI-210 Gas injector</td></tr> <tr><th>Sample</th><td>Non-pressurized gas, Volatile liquid</td></tr> <tr><th>Injection</th><td>10µl for liquid, 10ml for gas</td></tr> <tr><th>Carrier</th><td>Argon</td></tr> <tr><th>Heat</th><td>80°C for liquid</td></tr> <tr><th>Power</th><td>AC 100V to 240V, 50/60Hz, 20VA</td></tr> <tr><th>Dimensions / Weight</th><td>220(W) × 200(D) × 110(H) mm / 4 kg</td></tr> </thead> </table>	MODEL	GI-210 Gas injector	Sample	Non-pressurized gas, Volatile liquid	Injection	10µl for liquid, 10ml for gas	Carrier	Argon	Heat	80°C for liquid	Power	AC 100V to 240V, 50/60Hz, 20VA	Dimensions / Weight	220(W) × 200(D) × 110(H) mm / 4 kg		
MODEL	ASC-250L Liquid sample changer																																
Sample	Liquid (non-aqueous, aqueous)																																
Injection	Max 150µl (depend on sample)																																
Injection speed	0.4 - 50µl/sec (depend on sample)																																
Number	50pos in each 2, 4, 6ml vial tray																																
Power	AC 100V to 240V, 50/60Hz, 180VA																																
Dimensions / Weight	460(W) × 320(D) × 470(H) mm / 16 kg																																
MODEL	GI-210 Gas injector																																
Sample	Non-pressurized gas, Volatile liquid																																
Injection	10µl for liquid, 10ml for gas																																
Carrier	Argon																																
Heat	80°C for liquid																																
Power	AC 100V to 240V, 50/60Hz, 20VA																																
Dimensions / Weight	220(W) × 200(D) × 110(H) mm / 4 kg																																